

SPEC 4 DISINFECTANT WIPES #SP822802

Disinfects, Sanitizes, Deodorizes and Cleans



Features:

- Disinfects in just 4 minutes
- Kills 99.9% of bacteria in just 5 minutes
- Phenol free, alcohol free, bleach free and non-staining
- 180 wipes per canister - wipe size 6"x7"

Description:

SPEC4 WIPES utilize a powerful formula that is phenol free, alcohol free, bleach free and non-staining. Ideal for cleaning, disinfecting and sanitizing hard, nonporous, inanimate surfaces in many applications. Effective against a broad range of bacteria. Disinfects in only 4 minutes and sanitizes in only 15 seconds.



Active Ingredients:

Octyl decyl dimethyl ammonium chloride	0.0909%
Dioctyl dimethyl ammonium chloride	0.0364%
Didecyl dimethyl ammonium chloride	0.0545%
Alkyl (C14, 50%; C12, 40%; C16, 10%) dimethyl benzyl ammonium chloride.....	0.1212%

Effective Against:

- Influenza A/Brazil influenza
- Escheria coli
- Staphylococcus aureus
- Listeria monocytogenes
- Pseudomonas aeruginosa
- Herpes simplex I & II
- Klebsiella pneumonia
- Streptococcus pyogenes
- Hepatitis B & C
- Campylobacter jejuni
- Salmonella schottmuelleri
- Many others

For Use By:

- Hospitals
- Dental Offices
- Nursing Homes
- Correctional Facilities
- Daycare Centers
- Exercise Facilities
- Bathrooms
- Schools
- Factories
- Buses
- Terminals
- Hotels

Product Characteristics:

Appearance:	White towel
Odor:	Lemon
Vapor Pressure:	No data available
Specific Gravity:	1 g/ml using the liquid component of the towelette
Water Solubility:	Liquid component is soluble in water
Relative Density:	No data available
Flash Point:	>200.01°F Closed Cup using the liquid component of the towelette
VOC (%):	<1% using the liquid component of the towelette
pH:	10.5 – 12 (Tested using the liquid component of the towelette)
Storage:	Keep container tightly closed when not in use. Store in a cool area.
Transport Information:	
UN Number:	Not regulated for transport
Proper Shipping Name:	Not regulated for transport
Class:	Not regulated for transport
Packing Group:	Not regulated for transport

Directions:

Complete directions on product label.

To Sanitize/Deodorize: Treated surfaces to remain wet for 15 seconds.

To Disinfect: Treated surface to remain visibly wet for 4 minutes for most bacteria. Consult product label for directions.

Pictograms:



Signal Word:

WARNING.

Personal Protective Equipment Required:

Use appropriate personal protective equipment when risk assessment indicates this is necessary.

DOT Placard:

Not regulated for transport.